

3D Organon XR User Manual **Ultrasound**

for Standalone & PC VR Headsets

Available for:

Standalone VR headset: Meta Quest 2 | Meta Quest 3/3s | Meta Quest Pro | Pico Neo 3 Pro | Pico 4E | Pico 4 Ultra E | HTC Vive Focus 3 | HTC XR Elite | Apple Vision Pro

PC VR Headsets: HTC Vive Pro/Pro 2/Pro Eye | HTC Vive Cosmos and Windows MR headsets




Table of Contents

Kindly continue reading, or click on a link to get directly to a specific section.

<u>1. Requirements</u>	3
<u>2. The Lobby</u>	4
Menu	
<u>3. How to exit the Ultrasound module</u>	5
<u>4. Settings</u>	6
Get Started	
<u>5. Ultrasound Probes</u>	7
<u>6. How to use the Ultrasound Probe</u>	7
<u>7. Adjusting the bed height</u>	8
<u>8. Adjusting the Model Position</u>	8
<u>9. Gender Selection</u>	9
<u>10. Ultrasound Settings</u>	10
<u>11. Saving your Probe positions</u>	11
<u>12. Labelling and Color Mapping</u>	12
<u>13. Drawing Tool</u>	13
Remote Delivery / Medverse	
<u>14. Join a session</u>	14
<u>15. Create a session</u>	15
<u>16. Inside the virtual room</u>	16
<u>17. Network Services Documentation</u>	17

1 Requirements

- For Hardware Requirements & Installation Instructions

PC VR headsets: click [here](#)

Standalone VR headsets: click [here](#)

- To log in, you must have an active **Subscription** which is provided to you with a License Key.

If you **already have a Subscription** then:

- Activate your 3D Organon account from the **Welcome to 3D Organon** email you received.
- Get your License key from the **Order confirmation | 3D Organon License key(s)** email you received.


If you **do not have a Subscription**:

- Contact our Customer Development Team at prosupport@3dorganon.com


2 The Lobby features

Select from the following modules:

- Ultrasound



3 How to exit the Ultrasound module



Click here to go back to the Lobby

Settings

4 Settings

Themes Blue (default), Black, Green, White

Users can download additional languages.

Options

You can also access our user manuals and FAQs for more information by clicking on the 'Help' button.

Logout

Press the 'Logout' button to shut down the app and remove your access credentials from autologin (recommended for public computers and libraries to avoid credentials from being stored in the device).

Quit

Press the 'Quit' button to turn off the app (your credentials are stored locally. You do not have to insert them again when you log in the same device). Internet access at login is still required.



Italiano



Once a language has been downloaded, you can uninstall it by clicking on the Bin icon.



5 Ultrasound Probes

- Use the **Trigger button** to press virtual buttons and move scrollbars. It is also used for Labelling and Color Mapping.

The **Grip button** is used to grab and move the Probe.


- There are three different Probes available:

1. Curved
2. Linear
3. Cardiac




6 How to use the Ultrasound Probe

- Turn your controller facing downwards, click and hold the **Grip button** to simulate the way the probe is held and moved in real life.




7 Adjusting the bed height

- To adjust the bed height, press the top or bottom side of the touchpad or thumbstick of the left controller.




8 Adjusting the Model Position

- To adjust the scanning position (Sagittal and Coronal), press the top or bottom side of the touchpad or thumbstick of the right controller.



9 Gender Selection

- By pressing the 'Gender' icon with the trigger button, you can switch between male and female models.





10 Ultrasound Settings

- Slide the scrollbar left or right with the trigger button, to change the Depth, Noise Opacity and Attenuation settings.
- Click on 'Reset' to bring them back to the default settings.




11 Saving your Probe positions

- By pressing the 'Probe Positions' icon with the trigger button, you can save your probe positions so that you can easily find it when you re-enter the ultrasound module.



You can erase a saved probe position by clicking on the 'Delete' icon.


Click on 'Clear All' if you need to delete all of your previously saved probe positions.



Click on 'Save' and then enter a probe position name.

12 Labelling and Color Mapping

- Select and Color Map an organ or structure in the 3D Model (Window 2) to view the corresponding color mapped organ in the Sonogram (Window 1), and vice versa.
- The name of the Organ is displayed in the title window at the top. For example, the image below depicts the Colon (large intestine).



13 Drawing Tool

- Click on the 'Screenshot' icon to make the window appear.
- Select a color from the color palette and choose the size of your brush to begin drawing.

It is possible to adjust the thickness and color of the line.

Clear everything by clicking 'Clear All'.

PC VR: To take a screenshot, click the 'Save Screenshot' icon. Then, navigate to the Pictures folder in File Explorer on your PC to locate it.

Close the window

Colon (large intestine)

Brush Size: 5

Clear All

Save Screenshot

Close

Gender

Screenshot

The image shows a VR interface for a 3D anatomical model of the colon. A control panel is overlaid on the left side of the screen, featuring a 'Brush Size' slider set to 5, a color palette with six options (grey, green, cyan, yellow, red, purple), and a 'Clear All' button. Below this panel are 'Save Screenshot' and 'Close' buttons. On the right, a 3D model of the colon is displayed with a red dashed circle drawn around a section. A 'Screenshot' icon (a camera) is highlighted with a red box. The background shows a female avatar lying on a table in a VR environment. The text 'Colon (large intestine)' is visible at the top of the interface.


Remote Delivery / Medverse

14 Join a session

- Owners of 3D Organon apps can join a multi-user training session, deployed by an owner of a Premium license or Institution. You can also join a multi-user session in Guest mode.


Click on the 'Network' icon.

Select from the list of the available rooms and press 'JOIN'.
You can follow the session as presented by the tutor and communicate using text chat functions.



Network Rooms				
Room name	Host name	Module	Clients	Access
Medis	James	Anatomy	1/11	Enter JOIN

Click on the 'Profile' icon if you want to change your Username and Avatar.



If a password is required to join a particular room, please contact your institution/instructor.



Remote Delivery / Medverse

15 Create a session

- Owners of the **Premium edition** are also able to create virtual rooms and host remote delivery training sessions. Other users with access to the internet can join a session from anywhere in the world.

Click on the 'Network' icon.

Then, to create a room, click on the '+' button.



Enter the desired Room name, number of supported Clients, room visibility (Private/Public), and an optional Password. Press the 'CREATE SESSION' button to initiate your session.

Create Room

Room name

Clients


Private

Password

CREATE SESSION

16 Inside the virtual room

Users can type messages and ask questions in the chat room.



To exit the room and return to the Lobby, press the 'Leave Room' button.



17 Network Services Documentation

By default, the multi-user function in 3D Organon XR uses three essential services. These services include 3D models, voice chat, and text chat synchronization across users and educators in a virtual anatomy training session.

Ensure you allow the following ports and protocols for full access to 3D Organon networking services:

TCP: 843, 943, 4520, 4530-4532, 4540, 9090-9092

UDP: 5055-5058

Please ensure you allow traffic to and from the following domains:

license.3dorganon.com

def.3dorganon.com


download.3dorganon.com

Chinese clients need to white-list the following domain on their local network:

license-cn.3dorganon.com

Please ensure these ports are open in your router/firewall AND the firewall in your operating system (e.g., Windows firewall).

In case you are experiencing problems with either of these services, please contact the IT helpdesk in your institution. They will know your network infrastructure requirements and have access to modify firewall settings.



For any questions or feedback, please contact our Local Distributor
or reach us at: support@3dorganon.com

Check out our [Knowledge Base](#) for more information about 3D Organon, at any time!

